



---

**MALDI-TOF Mass Spectrometry:  
A New Rapid ID Method in  
Clinical Microbiology**

**Patrick R. Murray, PhD  
WW Director, Scientific Affairs  
BD Diagnostic Systems**



# Outline

- MALDI-TOF is the most important innovation in organism identification since microorganisms were first grown in the laboratory.
- By the end of this decade, laboratories will routinely use MALDI-TOF for the routine identification of all bacteria, mycobacteria, yeasts and molds.




# Mass Spectrometry

- Mass spectrometer is composed of 3 functional units: ion source, mass analyzer, and detection device
- **Ion source** – ionize sample molecules and transfers into gas phase
  - MALDI considered “soft ionization” (intact proteins)
- **Mass analyzer** separates ions by mass-to-charge ( $m/z$ ) ratio, e.g. time of flight (TOF)
- **Detection device** used to record separated ions


# MALDI-TOF Ionization of Sample

## Matrix Assisted Laser Desorption/Ionization




# MALDI-TOF Ionization of Sample

Matrix Assisted **Laser** Desorption/Ionization




# MALDI-TOF Ionization of Sample

Matrix Assisted Laser **Desorption/Ionization**



# MALDI-TOF Ionization of Sample

**Acceleration  
electrode**




The positively charged proteins drift in the flight tube; smaller proteins drift faster

Detector




Drift region




# MALDI-TOF MS profile spectrum

## Mass Range: 2-20 kDa




# Unique Profiles for Bacteria, Yeasts and Molds



# MALDI-TOF Mass Spectrometry

## Pre-analytical Processing



- Single colonies can be either transferred directly to the target plate or extracted with ethanol or methanol.
- If sufficient number of organisms are present, broth cultures (e.g., blood culture broths) can be concentrated and processed.
- Cell wall structure is disrupted with a strong organic acid (e.g., formic, trifluoroacetic, or acetic acid).
- Acetonitrile is used for protein extraction.
- Extracts spotted on stainless steel plates coated with Teflon (or disposable cards) and overlaid with matrix.



# Target Preparation - Rapid

## Direct Smear Method

- Touch colony with transfer device, such as toothpick
- Transfer a small amount onto spot
- Air dry
- Cover with 1  $\mu$ L of matrix; air dry
- Analyze



# Target Preparation - Best

## Modification of Extraction Method (Haigh et al, J Clin Microbiol 2011)

- Touch colony with transfer device, such as toothpick
- Transfer a small amount onto target plate; air dry
- Add 1  $\mu$ L of absolute **formic acid** to colony; air dry
- Cover with 1  $\mu$ L of matrix; air dry
- Analyze

# Target Preparation – Mycobacteria

- Suspend a small loopful of growth in water and heat to 95°C for 30 min to kill mycobacteria
- Disperse organisms with micropestle; centrifuge; remove supernatant
- Resuspend pellet in 70% formic acid and silica beads (breaks open the organisms); vortex 10 min
- Add acetonitrile (extracts proteins), vortex 10 min, centrifuge
- Apply supernatant to MALDI target plate; add matrix; analyze
- **Less than one hour from detection to ID**


# Target Preparation – Filamentous Fungi

Lau et al, J Clin Microbiol 2013


- Suspend a small colony fragment in absolute ethanol and silica beads; vortex 15 min.
- Remove ethanol by centrifugation; resuspend pellet in 70% formic acid and vortex 5 min
- Add acetonitrile, vortex for 5 min; centrifuge
- Apply supernatant to MALDI target plate; add matrix; analyze
- **Less than 30 min from detection to ID**

# MALDI-TOF Identification Protocol for Processing Blood Culture Broths

- Blood cells, serum proteins, and hemoglobin must be removed
- Lyse blood cells and concentrate bacteria by centrifugation.
- Wash pellet twice with sterile water to remove serum proteins
- Wash with 70% ethanol to kill bacteria
- Continue processing as with colonies from agar plates



# MALDI-TOF Profiles: Colonies vs. Broth



Ferroni et al, JCM 48:1542-48, 2010

# MALDI-TOF Microbial Identifications

- More than 250 publications in 2012
- Identification accuracy demonstrated for bacteria, mycobacteria, yeasts, and molds
- Accuracy determined by preanalytical processing, number of organisms analyzed, taxonomy of isolate, and database
- Most non-identifications due to:
  - too few organisms in sample
  - isolate not in database




# Effect of MALDI on Antibiotic Treatment of Bacteremic Patients

- Vlek et al (PLoS ONE, 2012) assessed the impact of MALDI-TOF on treatment of bacteremic patients:
  - Direct ID of blood culture isolates reduced time to ID by 28.8 hours
  - The proportion of patients receiving appropriate antibiotic therapy increased by 11.3%.
  - The effect was directly related to guided empiric therapy because time to definitive AST results was not modified in this study.

# Impact of MALDI, Rapid AST, and Antibiotic Stewardship

Perez et al, Arch Pathol Lab Med, 2013



## Length of Stay and Cost Outcomes

Outcome	Preintervention Cohort (n=100)	Intervention Cohort (n=101)	P value
Hospital length of stay	11.9 days	9.3 days	0.01
ICU length of stay	7.3 days	6.3 days	0.05
Total hospital costs	\$45,709	\$26,162	0.009

# Comparison of MALDI vs. Standard Identification Tests

Tan et al, J Clin Microbiol 2012

- Large medical center laboratory compared identification of all bacteria and yeasts isolated during a 12 week period by MALDI-TOF (Bruker) and standard ID tests in routine use (biochemical, agglutination, GLC, sequencing).
- 952 isolates from 2,214 specimens were processed.
- MALDI protocol provided ID an average of 1.45 days earlier.
- The estimated annual cost of organism identification was reduced by \$102,424 or 57%.

# MALDI: Workflow Impact

- MALDI can replace all biochemical and morphological identifications for bacteria, fungi, mycobacteria and blood isolates
  - Decreased consumable supplies
  - Decreased QC tests
  - Decreased repeat testing
  - Decreased supplemental tests
- Anaerobes can be identified as easily as aerobic bacteria
  - Elimination of subcultures to demonstrate strict anaerobes
  - Elimination of preliminary and definitive biochemical tests
  - Reduction to time to ID from days to minutes
- Definitive identification of molds and mycobacteria the same day growth is detected vs. days to weeks with traditional approaches